

Database Programming Project

Team 3: Xufang Deng, Sean Higgins, Yukimura Otake

Laboratory Management System

- What is organization for whom this project is being developed?

- Research laboratories


- What's the purpose and scope of this project?

- Improving documentation, reporting, and communication, etc.


- Who are the users of the project?

- laboratory personnel: Principal Investigator, Scientists, Lab manager

Laboratory personnel relationship









Laboratory Management System



Laboratory Database in phpMyAdmin

phpMyAdmin



Recent Favorites

+

information_schema

-

Laboratory

+ New

+ Employee

+ Equipment

+ Experiment

+ Product

+ Project

+ Reagent

+ Uses_Equipment

+ Uses_Reagent

laboratory_backup

Server: localhost:8


Structure SQL


Filters










Containing the word:









	Table
<input type="checkbox"/>	Employee
<input type="checkbox"/>	Equipment
<input type="checkbox"/>	Experiment
<input type="checkbox"/>	post
<input type="checkbox"/>	Product
<input type="checkbox"/>	Project
<input type="checkbox"/>	Reagent
<input type="checkbox"/>	user
<input type="checkbox"/>	Uses_Equipment
<input type="checkbox"/>	Uses_Reagent
	10 tables





Table structure and Relation view of Experiment table


 **Table structure**

 **Relation view**







	#	Name	Type	Collation	Attributes	Null	Default	Comments	Extra	Action
<input type="checkbox"/>	1	experiment_ID 	char(4)	utf8_general_ci		No	None			 Change
<input type="checkbox"/>	2	project_ID 	char(8)	utf8_general_ci		No	None			 Change
<input type="checkbox"/>	3	employee_ID 	int(9)			No	None			 Change
<input type="checkbox"/>	4	experiment_Objective	varchar(100)	utf8_general_ci		No	None			 Change
<input type="checkbox"/>	5	date	date			No	None			 Change
<input type="checkbox"/>	6	results	text	utf8_general_ci		Yes	NULL			 Change

 ☐ Check all With selected:  Browse  Change  Drop  Primary  Unique  Index  Fulltext

 [Print](#)  [Propose table structure](#)  [Move columns](#)  [Normalize](#)

Add column(s) 

Indexes

Action	Keyname	Type	Unique	Packed	Column	Cardinality	Collation	Null	Comment
 Edit  Drop	PRIMARY	BTREE	Yes	No	experiment_ID	3	A	No	
 Edit  Drop	Employee ID	BTREE	No	No	employee_ID	3	A	No	
 Edit  Drop	project ID	BTREE	No	No	project_ID	3	A	No	

Tue Dec 10 2019 20:23:38 GMT-0600 (Central Standard Time)

Search..

1 1 Identifying oncogene myc 2019-01-01

Identifying oncogene myc

2 2 Construct a plasmid DNA expressing myc oncoprotein 2019-01-02

Construct a plasmid DNA expressing myc oncoprotein

3 3 Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

Creating an oncogene knockout cell line with CRISPR-Cas

Our Sidebar

Useful links.

Reagent Inventory

Equipments

Experiment

Create a new experiment

New Experiment

Experiment ID

Project ID

LU201501

Employee ID

123456789

Objective

Experiment date:

mm/dd/yyyy

Results

Create this experiment

Create a new experiment

You have added a new experiment!

1 1 Identifying oncogene myc 2019-01-01

Identifying oncogene myc

2 2 Construct a plasmid DNA expressing myc oncoprotein 2019-01-02

Construct a plasmid DNA expressing myc oncoprotein

3 3 Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

Creating an oncogene knockout cell line with CRISPR-Cas

4 4 test1 2019-01-01

test1

Click individual experiment to display the details

Identifying oncogene myc 2019-01-01

Update

Delete

Add product

Add equipment

Add reagent

Experiment results: Obtained the entire sequence of oncogene myc

Experiment product(s):

- Product Name: oncogene myc cDNA
- Product Name: Purified P53 protein

Experiment equipment:

- Equipment ID: BR1001
- Equipment ID: TF1002

Experiment reagent(s):

- Reagent ID: R0220
- Reagent ID: TF1222

Delete an experiment

The experiment has been deleted!

1 1 Identifying oncogene myc 2019-01-01

Identifying oncogene myc

2 2 Construct a plasmid DNA expressing myc oncoprotein 2019-01-02

Construct a plasmid DNA expressing myc oncoprotein

3 3 Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

Creating an oncogene knockout cell line with CRISPR-Cas

Add a new product

You have added a new product!

Identifying oncogene myc 2019-01-01

Update

Delete

Add product

Add equipment

Add reagent

Experiment results: Obtained the entire sequence of oncogene myc

Experiment product(s):

- Product Name: oncogene myc cDNA
- Product Name: Purified P53 protein
- Product Name: test Name

Experiment equipment:

- Equipment ID: BR1001
- Equipment ID: TF1002

Experiment reagent(s):

- Reagent ID: R0220
- Reagent ID: TF1222

Add a new equipment

You have added a new equipment!

Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

- Update
- Delete
- Add product
- Add equipment
- Add reagent

Experiment results: Yes, I created

Experiment product(s):

Experiment equipment:

- Equipment ID: BR1001

Experiment reagent(s):

Add a new reagent

You have added a new reagent!

Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

Update

Delete

Add product

Add equipment

Add reagent

Experiment results: Yes, I created

Experiment product(s):

Experiment equipment:

- Equipment ID: BR1001

Experiment reagent(s):

- Reagent ID: R0220
-

Additional useful links and dropdown menu on the homepage

[Experiment Report](#) [Home](#) [About](#) [Menu](#) [New Experiment](#) [Account](#) [Logout](#)

Tue Dec 10 2019 20:23:38 GMT-0600 (Central Standard Time)

1 1 Identifying oncogene myc 2019-01-01

Identifying oncogene myc

2 2 Construct a plasmid DNA expressing myc oncoprotein 2019-01-02

Construct a plasmid DNA expressing myc oncoprotein

3 3 Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

Creating an oncogene knockout cell line with CRISPR-Cas

Our Sidebar
Useful links.

Reagent Inventory

Equipments

Experiment

[Menu](#) [New Experiment](#) [Account](#) [Logout](#)

Experiment

Equipments

Reagent

Look up inventory information is helpful for Scientists to choose what to use

[Catalog_number: R0220](#) Supplier: NEB Quantity: 2 Expiration_date: 2021-10-01

BamH I

[Catalog_number: TF1222](#) Supplier: Fisher Quantity: 4 Expiration_date: 2022-12-01

Ni purification column

[Equipment_ID: BR1001](#)

centrifuge

[Equipment_ID: TF1002](#)

PCR cycler

Cool stuff!

[Experiment Report](#) [Home](#) [About](#) [Menu](#) [New Experiment](#) [Account](#) [Logout](#)

Tue Dec 10 2019 20:23:38 GMT-0600 (Central Standard Time)

1 1 Identifying oncogene myc 2019-01-01

Identifying oncogene myc

2 2 Construct a plasmid DNA expressing myc oncoprotein 2019-01-02

Construct a plasmid DNA expressing myc oncoprotein

3 3 Creating an oncogene knockout cell line with CRISPR-Cas 2019-02-01

Creating an oncogene knockout cell line with CRISPR-Cas

Our Sidebar

Useful links.

Reagent Inventory

Equipments

Experiment

Use search bar to search experiment

2 Construct a plasmid DNA expressing myc oncoprotein 2019-01-02

Construct a plasmid DNA expressing myc oncoprotein

Code for the search bar

```
experiment_home.html x
<script src="https://ajax.googleapis.com/ajax/libs/jquery/3.4.1/jquery.min.js"></script>
{% extends "layout_laboratory.html" %}
{% block content %}
<input class="form-control" id="myInput" type="text" placeholder="Search..">
<br>
    {% for row in outString %}
        <article class="media content-section" id="inputs">
            <div class="media-body">
                <div class="article-metadata">
                    <a class="mr-2" href="#">{{ row['experiment_ID'] }}</a>
                    <small class="text-muted">{{ row['experiment_ID'] }} {{ row['experiment_Objective'] }} {{row['date']}}</small>
                </div>
                <h2><a class="article-title" href="{{ url_for('expt', experiment_ID=row.experiment_ID) }}">{{ row.experiment_Objective }}</a></h2>
            </div>
        </article>
        <script>
            $(document).ready(function(){
                $("#myInput").on("keyup", function() {
                    var value = $(this).val().toLowerCase();
                    $("#inputs *").filter(function() {
                        $(this).toggle($(this).text().toLowerCase().indexOf(value) > -1)
                    });
                });
            });
        </script>
    {% endfor %}
{% endblock content %}
```

Summary

- This application allows users (Scientists) to create, update, and delete an experiment
- Add additional info (reagent, equipment, etc.) to an experiment
- Allow to search experiments

This application 1) facilitate and motivate scientists to document their work; 2) make it easier for them to report; 3) also helpful for lab management, etc.

Future work

- Improve the display of each experiment;
- Create a request page so that Scientists can inform Lab manager to order items, etc.
- Complete the application for Principal investigator (PI), so PI can track the progress of each project
- Complete the application for Lab manager, so s/he can track inventory and manage the use of equipment